# SafeGuide<sup>®</sup> Over the Guidewire Esophageal Dilatation System

*Instructions for use and cleaning.  
Read carefully prior to use.*

Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Authorized Representative in EU:  
ADVENA LTD  
Pure Offices, Plato Close  
Warwick, CV34 6WE  
United Kingdom

063-0006W Rev. I, 2018/02, ECO 1895  
Made in the U.S.A.

## Diversatek Healthcare

### Corporate Headquarters

102 E. Keefe Ave.  
Milwaukee, WI 53212 USA  
P 800.558.6408 | 414.265.7620

[sales@diversatekhc.com](mailto:sales@diversatekhc.com)

**[DiversatekHealthcare.com](http://DiversatekHealthcare.com)**

**SAFEGUIDE® OVER THE GUIDEWIRE DILATATION SYSTEM**

*Instructions for Use and Cleaning.*

*Read carefully prior to use.*

**DESCRIPTION**

The SAFEGUIDE® OVER THE GUIDEWIRE ESOPHAGEAL DILATATION SYSTEM includes:

- A set of flexible, tapered, polyvinylchloride (PVC) dilators that have a central, longitudinal lumen for over-the-guidewire placement.
- A marked, stainless steel Spring Tip Guidewire, 210cm long with a removable end cap.
- Cleaning accessories such as a reusable brush and disposable cleaning adaptors.
- A storage case.

ALL OF THE ABOVE COMPONENTS MAY BE PURCHASED SEPARATELY.

**INTENDED USE**

For dilatation of upper esophageal webs, lower esophageal rings, caustic strictures, peptic esophageal strictures, and temporary ease of esophageal carcinoma.

The SafeGuide® Over the Guidewire Dilators are designed to be used over a pre-positioned SafeGuide® Guidewire, product numbers 1214-02 and 1214-02D.

**INSPECTION – PRIOR TO USE**

**CAUTION:** Remove protective end cap from new guidewire outside of procedure room and discard.

Inspect new dilators and accessories for damage from shipping and handling.

Inspect older dilators for wear. Worn dilators should NOT be used. Signs of wear can include: enlarged distal tip lumen due to repeated abrasions over the guidewire, surface cracks and crazing, cuts, and extremely faded print.

There is not an expiration date for the dilators because fitness for use is not age dependent. Fitness for use is determined by inspection for wear.

Visually inspect the guidewire for signs of wear or degradation such as kinks, bends, or breaks. Worn guidewires should NOT be used as they may not perform as expected.

**CONTRAINDICATIONS**

Contraindications include those specific to upper GI endoscopy.

Contraindications to dilation include, but are not limited to: uncooperative patient; asymptomatic strictures; inability to advance the dilator through the strictured area; coagulopathy; known or suspected perforation; severe inflammation or scarring near the dilation site, recent myocardial infarction, active ulcer and severe cervical arthritis.

**POTENTIAL COMPLICATIONS**

Potential complications associated with upper gastrointestinal endoscopy and esophageal dilation include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

**INSTRUCTIONS FOR USE**

1. Inspect all dilators and accessories before use. See INSPECTION - PRIOR TO USE section.
2. Perform screening endoscopy to identify strictured area.
3. Introduce the guidewire, spring tip end first, through the accessory channel of endoscope. A lubricant will help the guidewire glide through smoothly.  
Advance the guidewire until the spring tip is endoscopically visualized well beyond tip of scope and beyond the strictured area.  
The guidewire can also be externally monitored using the marking bands. The guidewire is marked in 20cm increments to help determine the location of the spring tip. The markings are groups of bands, beginning distally with a group of two (2) bands 40 cm from the tip. Each band represents 20 cm of distance to the spring tip, ending in a grouping of seven (7) bands 140 cm from the tip.  
A simple formula to measure the distance from a specific group of bands to the spring tip is:  
(number of bands) x 20 cm = distance in centimeters to the spring tip

Measurements are taken from the most distal band in each group

**CAUTION:** When placing the guidewire, discontinue advancement of the wire if resistance is met.

4. When the guidewire is in position well beyond strictured area, slowly begin to withdraw endoscope in 5-10cm increments while simultaneously advancing guidewire in 5-10cm increments to ensure guidewire remains in position.  
**CAUTION:** Continuous fluoroscopic monitoring of the guidewire is essential in order to ensure it remains in proper position.
5. When endoscope is removed completely, fluoroscopically confirm that guidewire has not been displaced.  
**Note:** The guidewire position can be externally monitored using the marking bands (as described in 2 above) in reference to the dental arch, both before and during dilatation.
6. Generously lubricate dilator and advance it over pre-positioned guidewire to strictured area.  
**Note:** SafeGuide® Dilators have two (2) series of numbers printed on each dilator on opposing sides. One series measures the distance to the distal tip, identified as **FROM TIP (American)**. The other series measures the distance to the stricture, identified as **FROM MAX O.D. (European)**. These printed distance markings can be used as a guide to avoid over-extending into the gastric cavity.
7. Proceed with esophageal dilation. Determine the SafeGuide® dilator to be used by the size of the strictures, starting with a smaller French sized dilator and gradually increasing the French size. Hold the dilator at the blunt, proximal end and place the lumen at the distal tapered tip over the end of the guidewire. Slide the dilator along the guidewire while continuing to monitor the guidewire location. Carefully remove dilator after completed dilation.  
Note: Successive dilations can be accomplished by repeating steps 5-6.
8. Upon completion of esophageal dilation, remove dilator and guidewire from patient.
9. Clean and prepare dilator for future use following cleaning instructions below.

#### **CLEANING INSTRUCTIONS**

##### **CAUTION:**

1. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, excessive wear, bends or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.
2. **DO NOT USE** Betadine to clean the dilators as it may cause discoloration of the material.
3. SafeGuide® Dilators **CANNOT** be steam autoclaved. The high temperatures of steam autoclaving will destroy the dilator.
4. **DO NOT USE** alcohol for cleaning or disinfection.
5. Appropriate gloves, gowns, and other protective wear should be worn when cleaning and sterilizing soiled or used dilators.
6. It is the responsibility of the user to qualify any deviations from the recommended methods of cleaning.

SafeGuide® Dilators must be thoroughly cleaned prior to use.

##### **PRE-CLEANING (SafeGuide® Dilators and Spring Tip Guidewires):**

1. Immediately after use, soak dilator/ guidewire in a solution of warm water and enzymatic cleaning solution. **Note:** Refer to the enzymatic cleaning solution instructions for recommended duration of cleaning.
2. Remove gross contaminates on the external surfaces of the dilator/guidewire by hand, using a soft bristle brush, sponge, or wipe, with warm water and enzymatic solution.
3. Clean the inner lumen of dilator using a solution of warm water and enzymatic solution and the SafeGuide® Cleaning Brush (either part number 1214-03D or 1214-03R). Insert the brush end first into the distal end of the dilator. Push the brush through the lumen, brushing the entire length of the dilator lumen.

4. Use the SafeGuide® Cleaning Adapters to flush the inner lumen with a solution of warm water and enzymatic cleaner.
5. Thoroughly rinse the dilator/guidewire (1214-02 only) with clean water. Use the SafeGuide® Cleaning Adapters to flush the inner lumen of the dilator with clean water.
6. Pat dry with clean, non-contaminating wipes. Use medical grade (oil-free) forced air to dry the lumen of the dilator.
7. Proceed with high-level disinfection / sterilization of dilator/guidewire (1214-02 only).

**HIGH-LEVEL DISINFECTION / STERILIZATION:**

The following disinfection solutions can be used with the Diversatek Healthcare SafeGuide® Dilators and accessories as indicated in the following table (marked with '✓'):

Solution	SafeGuide® Dilators	SafeGuide® Reusable Spring Tip Guidewire 1214-02	SafeGuide® Reusable Cleaning Brush 1214-03R
Glutaraldehyde	✓	✓	✓
Ortho-phthalaldehyde (OPA)	✓		✓
Peracetic acid	✓		
Hydrogen peroxide	✓		✓

Please refer to the manufacturer’s instructions for complete information and precautions for disinfection solutions with any of the above active ingredients.

The following specific brands have been validated with SafeGuide® Dilators and accessories as indicated in the following table (marked with '✓'). Please refer to manufacturer’s instructions for complete information and precautions:

System	SafeGuide® Dilators	SafeGuide® Reusable Spring Tip Guidewire 1214-02	SafeGuide® Reusable Cleaning Brush 1214-03R
Cidex® OPA High Level Disinfectant	✓		
Steris Resert® XL HLD High Level Disinfectant	✓	✓	✓
Rapicide PA High Level Disinfectant	✓		
Reliance™ EPS Endoscope Processing System	✓	✓	✓
Medivators Advantage Plus System	✓		
Medivators DSD series of automated endoscope reprocessors	✓		
Steris SYSTEM 1E® Liquid Chemical Sterilant Processing System	✓		
Steris V-Pro® Sterilization System		✓	✓

**Note:** Proper rinsing with clean water and complete drying of the dilator is necessary with any disinfection solution.

1. Completely immerse dilator in a disinfectant. **Note:** Refer to disinfectant solution instructions for duration of immersion required for high level disinfection or sterilization. If difficulty is encountered with solution penetration into the guidewire lumen, attach suction to the tapered end of dilator or utilize a SafeGuide® Cleaning Adaptor.

2. Rinse dilator thoroughly with clean, running water. Rinse the guidewire lumen by attaching suction to tapered end of dilator or using a SafeGuide® Cleaning Adaptor to suction water through this inner guidewire lumen. **Note:** Proper rinsing with clean water and complete drying of dilator and guidewire lumen is necessary following use of these solutions.
3. Dry external portions for dilator completely with a lint-free cloth. Use medical grade (oil-free) forced air to dry guidewire lumen.

The Reusable SafeGuide® Guidewire (without the protective end cap, part number 1214-02 only) and Reusable SafeGuide® Cleaning Brush (part number 1214-03R only) may be exposed to steam autoclave, using a standard steam sterilization cycle of 132°C for 15 minutes. The guidewires and reusable cleaning brushes should be individually wrapped for the steam sterilization cycle. After steam sterilization, allow the guidewires/cleaning brushes to completely cool to room temperature (25°C). Do not immerse in cold liquids. Dry thoroughly.

The Disposable SafeGuide® Cleaning brush (part number 1214-03D) and Disposable SafeGuide® Guidewire (part number 1214-02D) are single-use only. Do not reuse, reprocess, or resterilize single-use products. The materials used in the manufacture of the device may not withstand repeated reprocessing. The device may not perform as intended by the manufacturer if it is reused.

SafeGuide® Dilators are made with flexible, medical-grade polyvinylchloride (PVC). The products are compatible with ethylene oxide (EO) gas sterilization.

Device to be sterilized should first be placed in appropriate packaging materials capable of maintaining sterility, such as Tyvek® peel pouches. Typical conditions for ethylene oxide sterilization are 500 mg/L ethylene oxide, 50-70% relative humidity (RH) and 120-130°F with exposure time dependent on hospital regulations and/or sterilizer instructions. EO cycle temperatures should not exceed 135°F or 80% relative humidity (RH). EO cycle should have adequate aeration time to meet the requirements of ISO 10993-7.

The user must validate cycle efficacy along with the corresponding packaging/sterile barrier materials.

**STORAGE**

The SafeGuide® Dilators should be out of direct light, and at room temperature. SafeGuide® Dilators should be stored flat whenever possible. The SafeGuide® Storage Case provides a protected environment.

**DESCRIPTION OF SYMBOLS USED ON LABELS**

Manufacturer

Date of manufacture (YYYY-MM-DD)

Use-by date (YYYY-MM-DD)

Lot number

Part number

Non-sterile

Do not re-use (Part # 1214-02D, 1214-03D, 1214-04)

Consult instructions for use

Caution, consult accompanying documents

Not made with natural rubber latex

Contains Phthalate

Diversatek™ Healthcare is a trademark of Diversatek Inc.  
SafeGuide® is a registered trademark of Diversatek Healthcare.

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Diversatek™ Healthcare is een handelsmerk van Diversatek Inc.