

# CleanFreak® Endoscope Cleaning Sponge



Smart Cleaning Sponges and Pre-Cleaning Kit in convenient designs that ensure quick and confident pre-cleaning.

At Diversatek Healthcare, we look for ways to make your work easier. Our economical line of Cleaning Sponge designs and our Pre-Cleaning Kit ensure fast and reliable pre-cleaning.

# CleanFreak® Endoscope Cleaning Sponges



The pre-cleaning kit holds an enzyme-imbibed sponge in a water basin.

## Dry Imbided Sponges. Just add water.

Infused with a water-activated enzymatic cleaner, there's no need to inventory additional solutions. There's no mess with Diversatek Healthcare sponges compared to wet sponge packs.

## Multi-enzyme detergent.

To remove bio-burden, our imbided sponges contain protease, amylase and lipase as compared to other dual-enzyme competitors. Our dry imbided sponges have a low-foam formula for endoscope reprocessing.

## Durable cleaning solutions that fit the way you work.

If you prefer to use a specific cleaner, choose non-imbided sponges in cylindrical and flat textured designs. Whichever way you prefer to work, our economical sponges will help you get the job done.

### CleanFreak® Endoscope Cleaning Sponges Single-use

Part No.	Description	Color	Qty/Case
1201-01	Endoscope Sponge—Textured	Blue	200
1201-02	Endoscope Sponge—Cylindrical	Yellow	50
1201-03	Dry Imbided Endoscope Sponge—Cylindrical	Yellow	50
1201-20	Pre-Cleaning Kit—Imbided sponge, water basin	Sponge—Yellow Basin—Purple	24



The contoured design fully surrounds the endoscope for thorough cleaning.



Textured sponges conform to any product shape.

**DiversatekHealthcare.com**

P 800.558.6408 or 414.265.7620

**Corporate Headquarters**  
102 East Keefe Avenue  
Milwaukee, WI 53212 USA

orders@diversatekhc.com  
sales@diversatekhc.com

**Technical Research & Training Center**  
9150 Commerce Center Circle, Suite 500  
Highlands Ranch, CO 80129 USA

Learn more from  
your personal product  
representative at  
**800-558-6408**