

RePlay™ Hemostasis Clip



Clean deployments with continuous 360° rotation.

With continuous 360° rotation capability and clean separation during deployment, positioning the RePlay™ Hemostasis Clip is easy and precise.

RePlay™ Hemostasis Clip

Featuring the Smart Handle.

To facilitate the way you work, a simple twist of the positioning barrel rotates the clip 360°. Your fingers never have to be removed from the handle. The full rotation design ensures positioning without the unwanted “helicopter” effect.

Superior clip design provides clean deployment.

The clip separates cleanly from the sheath during deployment eliminating the need for additional shaking or pulling. Recognizable tactile feedback to the user helps reduce the risk of accidental deployment. In addition, the clip enables unlimited opening and closing prior to deployment.


Clip with confidence.

The geometry of the clip optimizes grasping power for enhanced vascular closure. The unique elastic telescopic design achieves hemostasis by focusing closure at an optimal tissue depth to maximize effectiveness.



RePlay™ Hemostasis Clip

Single-use, sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1170-01	11 mm width	230	2.8	5
1170-02	16 mm width	230	2.8	5



Rotating the positioning barrel on the handle turns the positioning clip smoothly with no “helicopter” effect.



Open and close the clip unlimited times prior to deployment for exact positioning.



Recognizable tactile feedback reduces the risk of accidental deployment.



Grasping power is optimized for better vascular closure.



The elastic, telescopic structure and clip geometry maximize the effect.

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters
102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center
9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408