


EndoGrip™ Grasping Forceps



Position confidently with smooth, 360° rotation.

EndoGrip™ Grasping Forceps position precisely while the rat tooth, alligator jaw ensures a powerful initial grab and secure retrieval.

EndoGrip™ Grasping Forceps

Featuring the Smart Handle.

The EndoGrip™ Grasping Forceps is designed to reduce movements. Without letting go of the handle, simply turn the positioning barrel to rotate the forceps 360° in either direction to retrieve objects from the gastrointestinal tract.

Secure, grab and hold design.


The versatile rat tooth design and serrated jaw holds onto the smallest surfaces for a confident grip.

EndoGrip™ Grasping Forceps Single-use, sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/ Case
1177-01	Rat Tooth-Alligator Combo	230	2.8	5



The powerful rat-tooth tip and serrated jaw design enable a secure grip on a wide variety of objects and surfaces.



Full, 360° positioning is achieved with a smartly designed rotating barrel that enables operation without letting go of the handle.



EndoGrip™ Grasping Forceps are ideal for both foreign body and stent removal. Grasp small objects of varying surfaces confidently.

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters
102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center
9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408