

Lasso™ Polypectomy Snare



Full, single-hand, 360° rotation simplifies control.

Lasso™ Polypectomy Snares from Diversatek Healthcare provide a smart, user-friendly solution to accommodate both hot and cold polypectomy.

Lasso™ Polypectomy Snare

The Smart Handle means responsive, precise control.


Rotate the positioning barrel left or right—a full 360°—without removing your other hand from the finger loops. Lasso™ Polypectomy Snares make your work easier by not interfering with the active cord.

Optimized flexibility.


Lasso™ Polypectomy Snares are designed for both hot and cold indications. Available in multiple sizes and loop geometries, the medium-stiff wire accommodates a wide range of user preferences for optimal flexibility.

Lasso™ Polypectomy Snares

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1180-01	Oval 10 mm	230	2.8	10
1180-02	Oval 15 mm	230	2.8	10
1180-03	Oval 25 mm	230	2.8	10
1180-04	Oval 32 mm	230	2.8	10
1180-10	Hex 15 mm	230	2.8	10
1180-11	3-in-1 Hex 9/18/25 mm	230	2.8	10
4525	Universal Active Cord	—	—	1



3-in-1 snare (1180-11) is adjustable to 3 distinct diameters to accommodate multiple polyp sizes, reducing the total number of snares needed in a procedure.



Smart, responsive 360° rotation is the gold standard in rotatable snares.



Multiple shapes and sizes to accommodate both hot and cold polypectomy.



True-to-size geometry of the snare designed to realize full intended size and maintain shape integrity after one or more deployments.

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters
102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center
9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408