


# M-Catch® Retrieval Nets



**Capture foreign bodies and polyps confidently.**

M-Catch® Retrieval Nets function exactly as you would expect with a crisp, responsive mechanism that feels like an extension of your own hand.



# M-Catch® Retrieval Nets

## Durable construction for multiple deployments.

M-Catch® Retrieval Nets feature netting material that is interwoven with the wire loop. This construction provides a secure assembly designed to be deployed multiple times. For a sure grip, the netting material conforms to the shape of the object being retrieved.

## Crisp, accurate deployment for confident retrieval.

The handle, sheath and responsive netting loop were designed to work together to enable crisp deployment for easier positioning around the target object. Available in three lengths, M-Catch® is matched to the endoscopes you use.


### M-Catch® Retrieval Nets

Single-use, non-sterile

| Part No. | Description             | Length (cm) | Minimum Working Channel (mm) | Qty/Case |
|----------|-------------------------|-------------|------------------------------|----------|
| 1175-00  | Standard oval, 3 x 6 cm | 160         | 2.8                          | 5        |
| 1175-01  | Standard oval, 3 x 6 cm | 230         | 2.8                          | 5        |
| 1175-02  | Standard oval, 3 x 6 cm | 350         | 2.8                          | 5        |
| 1175-10  | Polyp oval, 3 x 6 cm    | 230         | 2.8                          | 5        |



PTFE tubing resists kinking and bending while enabling smooth insertion through the endoscope working channel.



M-Catch® netting is interwoven with the wire loop and is designed for multiple deployments. Polyp netting option provides enhanced visibility for polyp retrieval.



Smooth actuation with a crisp deployment enables easy positioning around the target object.

[DiversatekHealthcare.com](http://DiversatekHealthcare.com)

P 800.558.6408 or 414.265.7620

**Corporate Headquarters**  
102 East Keefe Avenue  
Milwaukee, WI 53212 USA

orders@diversatekhc.com  
sales@diversatekhc.com

**Technical Research & Training Center**  
9150 Commerce Center Circle, Suite 500  
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com  
clinicalsupport@diversatekhc.com  
clinicaleducation@diversatekhc.com

Learn more from  
your personal product  
representative at  
**800-558-6408**