


BullDog® Disposable Endoscope Valves and Accessories



**Simplify endoscope use and clean up
with disposable valves and accessories.**

Cleaning an endoscope properly after a procedure is a critical part of the process. Cleaning simple components can be time consuming and difficult, and also presents the risk of cross-contamination. Reduce that risk with economical, single-use solutions from Diversatek Healthcare.

BullDog® Disposable Endoscope Valves and Accessories

Single-use components save time and reduce risk.

Reusable endoscope valves are difficult to clean effectively and not worth the time or risk. BullDog® single-use, disposable air/water and suction valves emulate the action of traditional OEM valves without the clean-up.

Engineered for a perfect fit.

The biopsy valve seals effectively as an instrument is passed down the endoscope's accessory channel. To accommodate larger instruments, simply lift the cap to reveal a larger entry port. And for convenience, color coding ensures you'll select the right valve matched to the endoscope brand.

BullDog® Disposable Endoscope Valves and Accessories

Part No.	Description	Qty/Case
1203-00	Biopsy valve, Olympus/Fujinon compatible, non-sterile	100
1203-05	Biopsy valve, Olympus/Fujinon compatible, non-sterile, bulk packed	200
1203-20	Biopsy valve, Pentax compatible, non-sterile	100
1203-50	Irrigating adapter, non-sterile	50
1203-70	2 piece kit: Air/water valve, suction valve, Olympus compatible, sterile	50
1203-71	3 piece kit: Air/water valve, suction valve, biopsy valve, Olympus compatible, sterile	50



Air/water and suction valve buttons feature a concave design that feels better during use. No lubricants are used that can cloud the lens and obscure the clinician's view.



Open the cap to accommodate larger endoscopic instruments.



Slit opening seals help prevent leakage and backflow.



Use the Irrigating Adapter to flush and clear the lens.



2- and 3-piece kits provide flexible ordering.

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters
102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center
9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408