


M-Catch™ Retrieval Nets



Capture foreign bodies and polyps confidently.

M-Catch™ Retrieval Nets function exactly as you would expect with a crisp, responsive mechanism that feels like an extension of your own hand.



M-Catch™ Retrieval Nets

Durable construction for multiple deployments.

M-Catch™ Retrieval Nets feature netting material that is interwoven with the wire loop. This construction provides a secure assembly designed to be deployed multiple times. For a sure grip, the netting material conforms to the shape of the object being retrieved.

Crisp, accurate deployment for confident retrieval.

The handle, sheath and responsive netting loop were designed to work together to enable crisp deployment for easier positioning around the target object. Available in three lengths, M-Catch™ is matched to the endoscopes you use.


M-Catch™ Retrieval Nets

Single-use, non-sterile

Part No.	Description	Length (cm)	Minimum Working Channel (mm)	Qty/Case
1175-00	Standard oval, 3 x 6 cm	160	2.8	5
1175-01	Standard oval, 3 x 6 cm	230	2.8	5
1175-02	Standard oval, 3 x 6 cm	350	2.8	5
1175-10	Polyp oval, 3 x 6 cm	230	2.8	5



PTFE tubing resists kinking and bending while enabling smooth insertion through the endoscope working channel.



M-Catch™ netting is interwoven with the wire loop and is designed for multiple deployments. Polyp netting option provides enhanced visibility for polyp retrieval.



Smooth actuation with a crisp deployment enables easy positioning around the target object.

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters
102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center
9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408