


Esophageal Dilator Storage



Store dilators easily in a space-saving,
wall-mounted storage cabinet.

Keep esophageal dilators clean, straight and organized between procedures
with an easy-to-install secure storage system.

Esophageal Dilator Storage




Easy to install compact, low profile storage cabinet system.

The Esophageal Dilator Storage system stores up to 24 SafeGuide® Dilators or WeightRight® Bougies just about anywhere.

Clean and secure.

Minimize exposure to contaminants through a forced air circulation system utilizing replaceable HEPA filters. A 405 nm LED disinfection light provides antimicrobial protection during storage.



Protect against bacterial and fungal antimicrobials during storage with an integral 405 nm disinfecting light.

Esophageal Dilator Storage

Part No.	Description	Qty/ Case
1214-70	Esophageal Dilator Wall-Mount Storage Cabinet. Includes fan with HEPA filter and 405 nm LED disinfection lighting system.	1

Replacement Parts & Accessories

1214-73	HEPA Filter	1
1214-74	Storage Cabinet Cleaning Brush	10
1214-75	Soft Grip Gasket	1
1214-76	Fan	1
1214-77	Power Supply	1
1214-78	LED Light Strip, Single Side	1
1214-79	Swing Bar	1



Minimize exposure to particulates and air-borne contaminants with forced air circulation and replaceable HEPA filters.



Medical grade silicone gently secures dilators vertically which keeps lumens clean.

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters
102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center
9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408